

Table 2-4 Tests and Inspection Requirements(Continued)

Test Block Name	Test or “Inspection” to Perform	Expected Test Results	What to Record on Service Record																		
	<p>Linearity Test</p> <table border="1" data-bbox="446 443 987 785"> <thead> <tr> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Reduce the manometer pressure down to 150 mmHg.</td> </tr> <tr> <td>2</td> <td>Wait 10 sec for the measurement to stabilize.</td> </tr> <tr> <td>3</td> <td>Compare the manometer values with the displayed values.</td> </tr> <tr> <td>4</td> <td>Document the value displayed by the Monitor. If the difference is greater than 3 mmHg then calibrate the module.</td> </tr> </tbody> </table> <p>Valve Test</p> <table border="1" data-bbox="446 894 987 1045"> <thead> <tr> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Raise the pressure to 280 mmHg.</td> </tr> <tr> <td>2</td> <td>Press Stop on the module to open valves.</td> </tr> <tr> <td>3</td> <td>Wait 5 seconds then document the value.</td> </tr> </tbody> </table>	Step	Action	1	Reduce the manometer pressure down to 150 mmHg.	2	Wait 10 sec for the measurement to stabilize.	3	Compare the manometer values with the displayed values.	4	Document the value displayed by the Monitor. If the difference is greater than 3 mmHg then calibrate the module.	Step	Action	1	Raise the pressure to 280 mmHg.	2	Press Stop on the module to open valves.	3	Wait 5 seconds then document the value.	<p>Value displayed on Monitor = x3 If difference (<= 3 mmHg)</p> <p>Proceed to next test</p> <p>Value displayed on Monitor = x4 (< 10 mmHg)</p>	<p>PN:P/x1/ x2/x3/x4 or PN:F/x1/ x2/x3/x4 where P = Pass and F = Fail</p>
Step	Action																				
1	Reduce the manometer pressure down to 150 mmHg.																				
2	Wait 10 sec for the measurement to stabilize.																				
3	Compare the manometer values with the displayed values.																				
4	Document the value displayed by the Monitor. If the difference is greater than 3 mmHg then calibrate the module.																				
Step	Action																				
1	Raise the pressure to 280 mmHg.																				
2	Press Stop on the module to open valves.																				
3	Wait 5 seconds then document the value.																				
<p>Performance Test Sidestream CO₂</p>	<p>Flow Adjustment Procedure</p> <table border="1" data-bbox="446 1173 987 1843"> <thead> <tr> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Enter the CO₂ Sidestream Task Window by pressing Parameters followed by CO2 Sidestream.</td> </tr> <tr> <td>2</td> <td>Press Start Pump. The date and time of the last flow adjustment are displayed.</td> </tr> <tr> <td>3</td> <td>Start the Flow Adjustment procedure by connecting the flow meter and pressing Start Pump. After about 5 seconds, the flow in ml/min displays in the Task Window.</td> </tr> <tr> <td>4</td> <td>Use the Adjust Flow to adjust the value displayed in the Task Window to the value displayed by the external flow meter.</td> </tr> <tr> <td>5</td> <td>Write down the value displayed by the CMS (x1).</td> </tr> <tr> <td>6</td> <td>Press Confirm to store the readjusted flow rate value. After about 5 seconds, the message “Adjustment done” displays. The flow is automatically set to 100 ml/min.</td> </tr> <tr> <td>7</td> <td>If the displayed value is not 100 ml/min, repeat Steps 3 through 5.</td> </tr> </tbody> </table>	Step	Action	1	Enter the CO ₂ Sidestream Task Window by pressing Parameters followed by CO2 Sidestream .	2	Press Start Pump . The date and time of the last flow adjustment are displayed.	3	Start the Flow Adjustment procedure by connecting the flow meter and pressing Start Pump . After about 5 seconds, the flow in ml/min displays in the Task Window.	4	Use the Adjust Flow to adjust the value displayed in the Task Window to the value displayed by the external flow meter.	5	Write down the value displayed by the CMS (x1).	6	Press Confirm to store the readjusted flow rate value. After about 5 seconds, the message “Adjustment done” displays. The flow is automatically set to 100 ml/min.	7	If the displayed value is not 100 ml/min, repeat Steps 3 through 5.	<p>Value displayed on CMS = x1 (where x1 = 100 +/- 10 ml/min)</p> <p>Proceed to next test</p>			
Step	Action																				
1	Enter the CO ₂ Sidestream Task Window by pressing Parameters followed by CO2 Sidestream .																				
2	Press Start Pump . The date and time of the last flow adjustment are displayed.																				
3	Start the Flow Adjustment procedure by connecting the flow meter and pressing Start Pump . After about 5 seconds, the flow in ml/min displays in the Task Window.																				
4	Use the Adjust Flow to adjust the value displayed in the Task Window to the value displayed by the external flow meter.																				
5	Write down the value displayed by the CMS (x1).																				
6	Press Confirm to store the readjusted flow rate value. After about 5 seconds, the message “Adjustment done” displays. The flow is automatically set to 100 ml/min.																				
7	If the displayed value is not 100 ml/min, repeat Steps 3 through 5.																				

Testing & Maintaining the CMS

Test and Inspection Requirements

Table 2-4 Tests and Inspection Requirements(Continued)

Test Block Name	Test or “Inspection” to Perform	Expected Test Results	What to Record on Service Record
Performance Test Sidestream CO ₂ (contd.)	Barometric Pressure Adjustment Procedure		Difference = x2 (<= 4 mmHg) PSC:P/x1/x2 or PSC:F/x1/x2 where P = Pass and F = Fail
	Step	Action	
	1	Press Barometer Pressure to adjust the Barometric Pressure value. The Task Window displays the stored barometric pressure in mmHg.	
	2	If the displayed value is incorrect, use the Barometer Pressure key to adjust the value to atmospheric pressure.	
	3	Document the Difference (x2) between the actual atmospheric pressure and the value displayed by the CMS.	
	4	Press Confirm to store the displayed or adjusted value. After about 5 seconds, the message “Adjustment done” displays. The barometric pressure is then set to the value you entered.	
	5	Press Main Screen (or Standard Display) to return to the standard display in Service mode	
Performance Test ECG	<p>Step 1 Connect the patient simulator to the ECG Parameter Module using the Patient cable.</p> <p>Step 2 Configure the Patient simulator as follows: ECG sinus rhythm HR = 120 BPM (Amplitude 1 mV)</p> <p>Step 3 Check displayed ECG wave and HR value against the simulator configuration.</p>	HR = 120 +/- 2 BPM	These results do not have to be reported.
Performance Test Respiration	<p>Step 1 Connect the patient simulator to the ECG/Resp Module using the patient cable.</p> <p>Step 2 Configure the patient simulator as follows: Base impedance line 1500 Ohm Delta impedance 0.5 Ohm Respiration Rate 40/min</p> <p>Step 3 Check displayed respiration rate against the simulator configuration.</p>	RPM = 40 +/- 2 /min	These results do not have to be reported.
Performance Test Cardiac Output	<p>Step 1 Connect the patient simulator to the C.O. Module using the patient cable.</p> <p>Step 2 Configure the Patient simulator as follows: Injection temperature: 2 °C Computation Const: 0.542 (Edward's Catheter) Flow 5 l/min</p> <p>Step 3 Check displayed value against the simulator configuration.</p>	C.O. = 5 +/- 1 l/min.	These results do not have to be reported.