

Calibration and Maintenance Schedule

The Atlas Monitor must be serviced by authorized Welch Allyn personnel or agents at 6 month intervals. Maintenance requirements are specified for 6 month and 12 month service intervals. The monthly CO₂ Reset operation can be performed by the user.

Service Interval

Maintenance Requirements

Every 6 months:

CO₂ Calibration

Every 12 months:

BP calibration, CO₂ reset, battery voltage calibration, printer adjustment, and temperature calibration.
Complete functional test

Table 2-2. Tools Required for Service.

| Description | Part# | Company |
|---|----------------------------|------------------------|
| 100cc Test Volume | T-112189 | Welch Allyn |
| 500cc Test Volume | T-112854 | Welch Allyn |
| T-10 TORX screwdriver | XTD-10 | Xcelite brand |
| 7/16" deep socket for Temperature Port nut | | Generic 1/4" drive |
| Squeeze Bulb and Valve | 5088-01 | Welch Allyn |
| Calibrated Manometer (0-10PSIG) | Digimano 1000 | Netech |
| Bio-Tek (NIBP TESTER) | BP Pump | Bio-Tek |
| Pneumatic Tubing | 5200-19 | Welch Allyn |
| Pneumatic Tubing (coiled) | 5200-19M | Welch Allyn |
| "Y" Fitting "Optional" | 9586TPK4 | Welch Allyn |
| "T" Fitting (3) | 9858TPK4 | Welch Allyn |
| Nonin Patient Simulator | 8000S | Nonin |
| Nonin Cable | 5200-52 | Nonin |
| Calibrated Thermometer for 90F to 115F range | 1002-3FC | ERTCO |
| LG, Adult Cuff/Bag | 5200-02 | Welch Allyn |
| Nellcor Patient Simulator For atlas with "MP 204/205 SpO2 PCB ONLY " | SRC-2 | Nellcor |
| Nellcor Patient Simulator "MP 506 PCB ONLY " | SRC-MAX | Nellcor |
| Nellcor Sensor Cable "Purple Connectors" will work on all Atlas units with 204/205/506 SpO2 PCB | DEC-8 | Nellcor |
| Nellcor Sensor Cable "Gray Connectors" "For Atlas units with MP 204/205 SpO2 PCB only" | EC-8 | Nellcor |
| ECG Simulator with Impedance Respiration | 214B | DNI Nevada |
| ECG Patient Cable (5 lead AHA) | 6200-02 | Welch Allyn |
| ECG Patient Cable (5 lead IEC) | 6200-04 | Welch Allyn |
| ECG leads, 5 Lead (IEC) | 6200-08 | Welch Allyn |
| ECG leads, 5 Lead AHA | 6200-06 | Welch Allyn |
| Certified Gas "10% CO2,10% O2, balance N2" | 0304724SRBD | Scott Medical Products |
| ETCO2 Water Trap (package of 5) | 6200-20 | Welch Allyn |
| ETCO2 Scrubber | 6200-21 | Welch Allyn |
| ETCO2 Adult Nasal Sample Line | 6200-22 | Welch Allyn |
| Surface Sensor, Temperature | 6200-15 | Welch Allyn |
| Atlas Repair and Calibration Software | 620538 | Welch Allyn |
| Printer Paper - Case | 6200-40 | Welch Allyn |
| Digital Multimeter with 10mV accuracy on a 10V scale; 10A Range | | |
| Adjustable DC power supply 5A @ 7V | | |
| Battery cable assembly | 620174-1 | Welch Allyn |
| Atlas interface cable to PC | 6200-60 | Welch Allyn |
| PC with Windows 95 and above and with HyperTerminal serial port software | Part of Windows 95 & above | |
| 1/4" Mono Phono Jack | | |
| 1000 Ohm precision resistor 1% | | |
| 1200 Ohm precision resistor 1% | | |
| 1350 Ohm precision resistor 1% | | |
| 1540 Ohm precision resistor 1% | | |
| 1870 Ohm precision resistor 1% | | |

Setting Date and Time

NOTE: Check date and time before doing any calibration. Set date and time if incorrect. Set time and date as follows;

1. Turn Atlas on.
2. Press **CLOCK** button next to power on button to check date and time. Use the far right **SELECT** button to scroll. Highlight the date or time that needs changed.
3. Press the **SET** button to adjust date and time.
4. Press **CLOCK** button to exit.

Pangea Communication Protocol

The Pangea Communication Protocol allows interaction of the Atlas with the computer through the serial interface port. A prompt is emitted at the computer screen when the instrument powers up and is ready to accept commands. The prompt is **Pangea>**. A prompt is emitted after the completion of each pangea command. Pangea commands are case sensitive. They are in English only.

BP Calibration

Required material.

- | | |
|------------------------------------|---------------------------------|
| 1. 500cc vessel | 4. Calibrated digital manometer |
| 2. Squeeze bulb with one-way valve | 5. Tubing and T fittings |
| 3. PC with HyperTerminal | 6. Serial cable |

NOTE: To start HyperTerminal you must have Windows 7 95 or higher installed on your computer.

1. From the Windows main screen follow the following sequence:
Start ⇒ Programs ⇒ Accessories ⇒ Communication ⇒ HyperTerminal
2. HyperTerminal setting are:
9600 baud rate, 8 bit word, 1 stop bit
no parity, no flow control
ANSI character set

- When you open HyperTerminal you will see a screen similar to the example in Figure 2-1. You will be prompted for a connection description. Choose any name. After you type in a name click **OK**.




Figure 2-1. Example of a new HyperTerminal connection.

- The next window you will see will be the window as shown in Figure 2-2. Click on **Connect using** then click on **COM1**.



Figure 2-2. Choosing COM1 in HyperTerminal.

5. Set the port settings as shown in the example in Figure 2-3. Now click **OK**.




Figure 2-3. HyperTerminal Port Settings.

6. Connect the pressure meter, bulb, and 500cc vessel to BP port with “T” connectors as shown in photograph Figure 2-4.
7. Connect the Atlas to PC with serial cable.
8. Turn the Atlas on. Start HyperTerminal on PC. Press the <Enter> key and you should see a **Pangea>** prompt.

NOTE: Take no more than 3 minutes for the 50mmHg calibration nor more than 3 minutes for the 250mmHg calibration as the Atlas will automatically, as a safety feature, open the blood pressure valve. If this happens you will have to turn the Atlas off then back on again and restart the calibration again.




Figure 2-4. Photo of Atlas BP calibration setup.

50mmHg Calibration

1. Enter the following commands at the Pangea prompt.
Pangea> **bp valve close** <ENTER>
Pangea> **bp safety off** <ENTER>
Pangea> **bp cal 5000** *Do not press* <ENTER> *yet!*
2. Raise the pressure with bulb to as close to 50.00mmHg as possible. Now press <ENTER>.
3. Release the pressure.

250mmHg Calibration

1. Enter the following command
Pangea> **bp cal 25000** *Do not press* <ENTER> *yet!*
2. Raise the pressure with bulb as close to 250.00 mmHg as possible. Now press < ENTER >
3. Enter the following command to save the calibration in the Atlas.
Pangea> **nvrw write** <ENTER>
4. Release the pressure.

CO₂ Reset (623xx Models Only)

- Required material.**
1. Watertrap
 2. Scrubber

NOTE: *The Scrubber looks similar to a watertrap, but it is filled with white granules. The scrubber is included with the 623XX models only.*

NOTE: *Make sure date and time are correct before performing the CO₂ reset.*

1. Turn Atlas on. Make sure the watertrap and scrubber are **NOT** attached to the Atlas.
2. Press the **DATE/TIME** button on the lower right of the monitor. The Set Date and Time and Other Options menu will be displayed.
3. Press the **CO₂/RESP ALARMS Off** button. The CO₂ Reset screen will appear.
4. You will see the following messages on the CRT.
 “CO₂ Reset”
 “Install CO₂ scrubber”
 “Press Trend to abort”
5. Install the watertrap to the Atlas. Install the scrubber to the watertrap.
6. You will see the following instructions on the screen.
 “Warming up” will be flashing on CRT.
 “May take up to five minutes” on CRT.
 “Press Trend to abort” on CRT.
7. After about 5 minutes you will see on the CRT.
 “CO₂ Reset”
 “Reset complete”
 “Remove CO₂ scrubber”
 “Press the trend button to exit”
8. Remove the CO₂ watertrap and scrubber.
9. Press **TREND** button to return to idle screen.



Replace watertrap after every six hours of use. Treat watertrap and used CO₂ sample lines as bio hazard material!

ET CO₂ Calibration

Required material:

1. Tank of approximately 10% CO₂, balance N₂ (certified) Blood Gas Mixture.
 2. Tubing and T connectors.
 3. Watertrap and scrubber.
1. Make sure the watertrap and scrubber are not attached to the Atlas. Turn the Atlas on.
 2. Place the instrument into the Service Mode by pressing the **DATE/TIME** button. Make sure date and time are correct. Press the **LEAD SELECT** button.
 3. Press **SELECT** button and scroll down to Calibrate CO₂.
The message “**Install CO₂ Scrubber**” will appear on the right side of screen.
 4. Attach the scrubber to the water trap.
 5. Insert the scrubber/water trap assembly into water trap socket. The message “**Enter span gas value using Set button 10%**” will appear.
 6. Press the **SET** button to change the value of span gas being used. The factory default value is 10%. Calibrate with a 8% to 12% certified CO₂ concentration known to be ±0.01%).
 7. The message press “**BP Start/Cancel**” will appear at the bottom right of CRT. Press the **BP/Start/Cancel** button.
 8. If you receive a “**Calibration Failed**” message at this point, check the date. If date is 2022 or above it will fail CO₂ calibration.
 9. Next you will see a message “**Warming up**”. After the Atlas warms up you will see a message “**Attach CO₂ gas**”. Remove the scrubber from the CO₂ water trap. Do not remove water trap.
 10. Attach the certified source of CO₂ gas to the CO₂ side-stream sampling tube as per Figure 2-5 below.

Adjust Regulator to
Approximately 2 psi




Figure 2-5. Example of CO₂ calibration setup.



CAUTION: IMPROPER USE, STORAGE OR HANDLING OF COMPRESSED GAS VESSELS CAN CAUSE DEATH OR INJURY. FOLLOW GAS MANUFACTURER'S SAFETY PROCEDURES!

11. Adjust the CO₂ regulator just enough to allow a small amount of gas to flow out of the vent (approximately 2 psi).
12. Press the **BP START/CANCEL** button. The message: "Sampling" will appear on CRT. "CO₂ calibration successful" or "CO₂ calibration failed" will appear on the CRT display.
13. Press the **TREND** button to exit.

No-Load Battery Voltage Calibration

NOTE: No-Load Battery Voltage Calibration procedure is for models 622 and 623 only. If you disconnect the battery you must reset time and date after you reconnect the battery.

Specifications: No load battery charge 6.85 VDC

Required materials: DVM

1. Remove the battery and disconnect the battery leads from the Atlas. Use the DVM to measure across the connectors, red+ and black-.
2. Adjust the voltage to 6.85 VDC by turning potentiometer R338, located behind the right battery jack. Turning clockwise will increase the no-load voltage and counter clockwise will decrease the no-load voltage.

Battery Voltage Calibration

Required materials.

1. DC power supply rated: 7 VDC at 5A
2. DMM / DVM with 10mV resolution on a 10 DC volt scale.
3. PC with HyperTerminal
4. Serial interface cable

1. Connect the serial cable to the PC and Atlas.
2. Hyper Terminal Settings are:
 - 9600 Baud, 8 bit word, 1 stop bit
 - no parity, no flow control
 - ANSI character set
 - Find HyperTerminal in Windows 95 or higher
 - Start ⇒ Programs ⇒ Accessories ⇒ Communication ⇒ HyperTerminal

NOTE: Make sure the Atlas **IS NOT** plugged into AC for this calibration procedure.

3. Remove and disconnect the battery from Atlas.
4. Set the power supply to 6.8VDC ± 200mV and connect the power supply to the battery connector on the Atlas.
5. Turn Atlas on..
6. Reduce the power supply to 6.0VDC and measure the voltage at the battery connector (at the Atlas) to the nearest 10mV.

NOTE: Do not measure at the power supply, since cable resistance will introduce error.

7. At the HyperTerminal prompt type:
Pangea> **power cal XXXX** <ENTER>

NOTE: XXXX represents the measured voltage in millivolts no decimal point. For example, if you measured 6.010VDC at the battery connector, use the command “**power cal 6010** <ENTER>”.

8. The Atlas will respond: raw = ZZZZ mV true = 6010 mV OK

NOTE: ZZZZ is the raw uncalibrated reading that the instrument made.

9. Reduce the power supply to 5.6 VDC. You should soon hear the “**low battery**” alarm.
10. Measure the voltage at the battery connector to the nearest 10mV.

11. At the HyperTerminal prompt type:

Pangea> **power cal XXXX** <ENTER>

NOTE: *XXXX represents the measured voltage in millivolts with no decimal point. For example if you measured 5.590 volts at the battery connector, then you would enter the command “**power cal 5590**” <ENTER>.*

12. The Atlas will respond:

raw = ZZZZmV true = 5590 OK

NOTE: *ZZZZ is the raw uncalibrated reading that the instrument made.*

13. Finish the calibration by typing:

Pangea> **hw reset** <ENTER>

NOTE: *This will re-boot the Atlas.*

14. Turn the Atlas off and reinstall the battery.

15. Reset time and date.

Printer Print Adjustment

1. Install new paper.
2. Turn Atlas on.
3. Attach an ECG simulator to Atlas and set simulator for a heart rate of 60 bpm, normal sinus rhythm.
4. Press **PRINT** button. Evaluate the darkness of waveform and text printout.
5. If either need to be changed press **DATE/TIME** button then press **LEAD SELECT** button to access Advanced Configuration menu.
6. Press **SELECT** button and scroll down to Printer Test Pattern. Then Press **HR ALARMS OFF** button
7. Two lines will be displayed:
 1. Waveform+128
 2. Text+ 78

NOTE: *These two numbers are the factory defaults and are a good starting point if the system is printing poorly or not at all.*

- The left **SET** button controls the waveform darkness and the right **SET** button controls the text darkness.

NOTE: Pressing the **SET** button up will increase the number and darken the waveform while pressing the **SET** button down will decrease the number and will lighten the waveform.

- Make changes to the printout as needed.
- Press **Trend** button to return to waveform screen.

Temperature Calibration

Required Material:

- | | |
|--|------------------------------|
| 1. PC with Windows7 95 or higher | 4. 1/4" mono phono jack |
| 2. Atlas serial cable | 5. Soldering iron and solder |
| 3. 1k Ohm, 1/2 watt precision resistor | 6. Ohm Meter |
- Solder the 1 K Ohm resistor to 1/4" mono phono jack.
 - Measure the resistance at the tip of the phono jack. Record that resistance reading, to two decimal points.
 - Plug the phono jack into the Atlas.

NOTE: The Atlas will show a temperature reading in the temperature display.

- Start HyperTerminal on PC.
- Hyper Terminal Settings are:
 - 9600 Baud, 8 bit word, 1 stop bit
 - no parity, no flow control
 - ANSI character set
 - Find HyperTerminal in Windows7 95 or higher
 - Start ⇨ Programs ⇨ Accessories ⇨ Communication ⇨ HyperTerminal
- At the Pangea prompt Type: PANGEA> **temp cal XXXXXX** <ENTER>

NOTE: XXXXXX is the resistance reading you measured and recorded in milliohms at the tip of the mono phono jack.

Example: If you measure 1000.40 ohms at the tip of the phono jack then you would type **temp cal 100040** <Enter>.

- Wait four seconds then type: PANGEA> **temp state** <ENTER>. You will see a value returned at the Pangea prompt.

8. Verify that the resistance given by the above command returns a value ± 0.5 ohms.
9. Verify temperature accuracy as outlined in Chapter 3, Table 3-2.

Calibration Date Set

NOTE: After calibrating the Atlas you must reset the calibration date. The calibration date is the date you performed the calibration. The calibration date appears in the Service Mode menu. See Figure 2-6 for an explanation of the Service screen.

1. At the Pangea prompt type.

```
Pangea> nvram set cal_date XXXX <ENTER>
```

NOTE: XXXX is the number of days from January 1, 1998 until the present date. See APPENDIX E for that number or calculate manually.

Explanation of an Atlas Service Screen




Figure 2-6. Example and Explanation of an Atlas Service Screen

Line A Indicates that the Atlas is a model 220 and that the Atlas has a printer. If the Atlas did not have a printer then the text “with printer” would be absent. The date field following indicates the date the Atlas was last calibrated.

Line B starts with a numeric sequence. The first three digits are the model number of the Atlas. The next five numbers are the Atlas serial number. The next sequence, alphanumeric, indicates what software version is currently loaded into the Atlas. The third sequence, alphanumeric, indicates what version boot software is currently loaded into the Atlas.

Line C starts with the SpO₂ OEM board used in the Atlas. There are two SpO₂ OEM boards used in the Atlas. One is Nellcor and the other is Nonin. The next sequence indicates the model of the SpO₂ board. The next sequence, that starts with a letter “V”, is the version software used with the current SpO₂ board. The date following the SpO₂ software is the date the OEM loaded the software into the SpO₂ board. SpO₂ OEM software can not be upgraded. If the most current software is needed, you will need to replace the SpO₂ board.

Line D is the CO₂ information if your Atlas has CO₂ installed. **Only models 623 have CO₂ installed.** If your Atlas does not have CO₂ installed then line **D** will be absent. If your Atlas has CO₂ installed then line **D** will start with CO₂ followed by a numeric number. That numeric value is the software loaded into the CO₂ board. The next alphanumeric sequence, starting with a “V”, is the version of that software. The next numeric sequence starting with a “#” is the serial number of the CO₂ board. The next sequence, a date, is the date the CO₂ was last calibrated. The last sequence, a date, is the date the CO₂ was reset. CO₂ OEM software can not be upgraded. If the current software is needed, you will have to replace the CO₂ board.

Table 2-3. Software Revision Table.

| Model | Operating System | Boot loader | Nellcor | Nonin | Pryon |
|-------|----------------------|---------------------|--------------------|-------|-----------|
| 621S0 | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | | V7 | |
| 621S0 | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | | V7 | |
| 621S0 | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | | V7 | |
| 621S0 | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | | V7 | |
| | | | | | |
| 621N0 | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | V1.8.1.0, 10/14/02 | | |
| | | | | | |
| 621SP | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | | V7 | |
| 621SP | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | | V7 | |
| 621SP | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | | V7 | |
| 621SP | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | | V7 | |
| | | | | | |
| 621NP | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/02 | V1.8.1.0, 10/14/02 | | |
| | | | | | |
| 622S0 | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | | V7 | |
| 622S0 | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | | V7 | |
| 622S0 | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | | V7 | |
| 622S0 | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | | V7 | |
| 622S0 | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | | V7 | |
| | | | | | |
| 622SP | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | | V7 | |
| 622SP | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | | V7 | |
| 622SP | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | | V7 | |
| 622SP | BB.02.5000, 6/14/02 | AA.01.0000, 6/20/99 | | V7 | |
| | | | | | |
| 622N0 | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | |
| 622N0 | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | |
| 622N0 | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | |
| 622N0 | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | V 1.8.1.0 10/14/02 | | |
| | | | | | |
| 622NP | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | |
| 622NP | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | |
| 622NP | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | V 1.2.0.0 10/14/02 | | |
| 622NP | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | V 1.8.1.0 10/14/02 | | |
| | | | | | |
| 623SP | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | | V7 | 0.E V1.00 |
| 623SP | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | | V7 | 0.E V1.00 |
| 623SP | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | | V7 | |
| 622SP | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | | V7 | |
| | | | | | |
| 623NP | AA.01.4000, 9/8/99 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | 0.E V1.00 |
| 623NP | BB.02.0000, 4/4/01 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | 0.E V1.00 |
| 623NP | BB.02.2000, 6/4/02 | AA.01.0000, 6/20/99 | V 1.2.0.0 12/17/97 | | |
| 623NP | BB.02.5000, 10/14/02 | AA.01.0000, 6/20/99 | V 1.8.1.0 10/14/02 | | |

Software Upgrade Procedure

NOTE: The following procedures are required to upgrade the software on a fully functioning Atlas or to reload software after replacing the CPU board. The download utility “**atlas_dl.exe**” loads the following files automatically.

1. atlas.out.gz
2. nvram_common.txt
3. nvram_(model#).txt
4. nvram_(language).txt

“Atlas_dl.exe” will also query the Atlas to determine what model number the Atlas is and what language to download.

Equipment or supplies required:

1. PC with Windows 7 95 or higher
2. Atlas serial cable
3. File: atlas_dl.exe (Included in the Atlas Repair Software listed in Table 2-2.)

NOTE: Make sure you have HyperTerminal turned off or the following utility download will not work!

1. Run the program **atlas_dl.exe** from the CD or copy the file to your hard drive and run the program from there.
2. Connect the serial cable between the Atlas and the PC’s “COM1 port”.
3. Turn your computer on and locate the “atlas_dl.exe” file.
4. Double left click on the “atlas_dl.exe” file.
5. When the file starts to download, the CRT will go blank on the Atlas.

NOTE: Do not use the computer while the program is downloading.

6. After downloading is complete, check all alarm settings and all user advanced configuration settings since these are **RESET** by this utility software download procedure.

NOTE: **Stop here** if you are just upgrading software on a fully functional Atlas monitor or if you have replaced the CPU board.

NOTE: If you have replaced the **MAIN BOARD** then you **MUST** continue with the next procedure (**DOWNLOADING NVRAM FILES**).

NOTE: Perform the following NVRAM downloading procedure if you have replaced the **MAIN BOARD**. The NVRAM resides on Main Board.

Down Loading NVRAM Files

Required Materials.

1. Computer with Windows 7 95 or higher.
 2. The latest Atlas repair software. Call your local Welch Allyn representative.
1. Once you have procured the latest Atlas Repair Software, you will be able to run the programs straight from the CD or copy the files to your hard drive.

NOTE: Through HyperTerminal you will need to load the following files:

1. *nvrाम_cal_init.txt*
2. *nvrाम_common.txt*
3. *model# of Atlas*
 - 3.1. You will load "*nvrाम200.txt*" if your Atlas is a model 621xx.
 - 3.2. You will load "*nvrाम210.txt*" if your Atlas is a model 622xx.
 - 3.3 You will load "*nvrाम220.txt*" if your Atlas is a model 623xx.
4. "*language.txt*" Any combination of or all of the following files
 - 4.1 "*nvrाम_english.txt*" for the English language
 - 4.2 "*nvrाम_french.txt*" for the French language
 - 4.3 "*nvrाम_german.txt*" for the German language
 - 4.4 "*nvrाम_spanish.txt*" for the Spanish language
 - 4.5 "*nvrाम_potuguese.txt*" for the Portuguese language
 - 4.6 "*nvrाम_italian.txt*" for the Italian language
 - 4.7 "*nvrाम_chinese.txt*" for the Chinese language
 - 4.8 "*nvrाम_japanese.txt*" for the Japanese language
5. *printer.txt* (if fitted with a printer), or *no_printer.txt*. (if not fitted with a printer)

NOTE: After loading the preceding files you, must then type the following commands at the Pangea prompt.

1. "*nvrाम set serial xxxxx*" where xxxxx is the serial number of the Atlas monitor.
 2. "*nvrाम write*" writes the information to memory.
 3. "*hw reset*" performs a hardware reset.
2. Connect the serial cable between the Atlas and the PC COM1 port.
 3. Open HyperTerminal program on PC.
Hyper Terminal Settings are:
9600 Baud, 8 bit word, 1 stop bit
no parity, no flow control
ANSI character set
Find HyperTerminal in Windows 7 95 or higher
Start ⇒ Programs ⇒ Accessories ⇒ Communication ⇒ HyperTerminal

4. Turn Atlas on. You should see the Pangea prompt. See example in Figure 2-7.




Figure 2-7. Pangea prompt in HyperTerminal.

5. Scroll over to **“Transfer”** and then scroll down and choose **“Send Text File”**. See example in Figure 2-8.




Figure 2-8. Choosing “Send Text File” in HyperTerminal.

6. Another window will then appear and prompt you for the location of the files. Double left click on that folder to open that directory. You can also run these programs from the CD. See example in Figure 2-9.




Figure 2-9. Choosing directory where NVRAM files are stored.

7. Once the folder is open, open the file “**nvr_{am}_cal_init.txt**” file by double left clicking on that file or by high lighting the file and then click on the **Open** button. See example in Figure 2-10.



Remember that after you have download the files in this section and you must complete this entire “Downloading NVRAM files**” section, that you must perform a complete calibration on the Atlas monitor.**




Figure 2-10. Choosing the nvr_{am}_cal_init.txt file.

NOTE: After the *nvram_cal_init.txt* has executed you should see a Pangea screen similar to window as shown in Figure 2-11.




Figure 2-11. Pangea screen after opening the *nvram_cal_init.txt* file.

8. From your pangea window choose **Transfer** then choose **Send Text File**.
9. Open the directory where the file *nvram_common.txt* is located or open from the CD.
10. Double left click on *nvram_comm.txt* file or highlight the file and then choose **Open**.

NOTE: After opening the *nvram_common.txt* file you should see a screen similar to the window shown in Figure 2-12.




Figure 2-12. Pangea screen after opening the *nvram_common.txt* file.

NOTE: Next you will download the language(s) that you want the Atlas to store in its Advanced Configuration menu for languages. You can have one or all the languages loaded in the Atlas. Listed in Table 2-4 are the language(s) choices and the files you will need to download to have the language(s) loaded in the Atlas. Each language you want loaded in the Atlas will require that the file associated with that language be loaded.

Table 2-4. Files needed to download language(s).

| Language | File needed to download |
|------------|-------------------------|
| English | nvrnm_english.txt |
| French | nvrnm_french.txt |
| Spanish | nvrnm_spanish.txt |
| German | nvrnm_german.txt |
| Portuguese | nvrnm_portuguese.txt |
| Italian | nvrnm_italian.txt |
| Chinese | nvrnm_chinese.txt |
| Japanese | nvrnm_japanese.txt |

13. Open the directory where the language files are located or open from CD.

14. Double left click on language.txt file or highlight the file and then choose **Open**.

NOTE: After you download a language.txt file you should see a Pangea screen similar to Figure 2-14.

```

Atlas - HyperTerminal
File Edit View Call Transfer Help
# Pangea> # Common NVRAM Configuration Settings for Models 200, 210 and 220.
Pangea> #
Pangea> noecho
Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea>
Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea>
Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea>
Pangea> # End of file
Pangea> -
Connected 0:00:25 Auto detect 9600 8-N-1 SCROLL CAPS NUM Capture

```

Figure 2-14. Pangea screen after loading a language file.

NOTE: Next you tell the Atlas if it does or does not have a printer.

1. If the Atlas has a printer you will download the file **nvr_{am}_printer.txt**.
2. If the Atlas does not have a printer then you will download the file **nvr_{am}_no_printer.txt**.

15. Open the directory where the printer.txt files are located or open from CD.

16. Double left click on the printer file or the no_printer file or highlight the printer file or the no_printer file and then choose **Open**.

NOTE: After you have downloaded the printer or no_printer file you should see a Pangea screen similar to the window as shown in Figure 2-15.

```

Atlas - HyperTerminal
File Edit View Call Transfer Help
#
Pangea> # NVRAM default printer settings for instruments WITH PRINTERS.
Pangea> #
Pangea>
Pangea> noecho
Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea>
Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea> Pangea>
Pangea> Pangea> Pangea>
Pangea> # End of file
Pangea>
Connected 0:00:18 Auto detect 9600 8-N-1 SCROLL CAPS NUM Capture

```

Figure 2-15. Pangea screen after downloading the nvr_{am}_printer.txt file.

NOTE: The next three commands will require that you actually type the command at the Pangea prompt.

NOTE: XXXXX on line 17 is the Atlas 5 digit serial number. The serial number is located on the bottom of the Atlas. Make sure that you have a space between the command **serial** and the five numbers you are entering.

17. At the Pangea prompt type the following command “nvr_{am} set serial XXXXX”<ENTER>

NOTE: After you have typed the **nvr_{am} set serial XXXXX** command and hit **ENTER** you should see a Pangea screen similar to the example in Figure 2-16.




Figure 2-16. Pangea screen after you have downloaded the nvrAm set serial command.

NOTE: To finish downloading the NVRAM files, type the next two commands at the Pangea prompt.

18. At the Pangea prompt type the command **nvrAm write** <ENTER>.

19. At the Pangea prompt type the command **hw reset** <ENTER>.

NOTE: The **hw reset** command will reboot the Atlas.



Remember that after you have download the files in this section and you **must** complete this entire “**Downloading NVRAM files**” section, that you **must** perform a complete calibration on the Atlas monitor.

NOTE: The firmware is not up grade able on the OEM boards (SpO₂ and CO₂) boards. If a higher version software is needed, then you will need to replace the OEM board.