| Monitor<br>Software | FMS-8 (M8048A), FMS-4 (865243, starting with H.15) |     |     |     |     |     |     |     |     |     |     |     |
|---------------------|--|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
|                     | B.0  | B.1 | C.0 | D.0 | E.0 | F.0 | G.0 | H.x | J.0 | K.2 | L.0 | M.0 |
| B.01)               | 0  | 0   | 0   | 0   | 0   | ×   | ×   | ×   | ×   | *   | ×   | ×   |
| B.11)               | 0  | 0   | 0   | 0   | 0   | *   | ×   | *   | ×   | *   | ×   | ×   |
| C.0                 | 0  | 0   | 0   | 0   | 0   | ×   | ×   | *   | ×   | *   | ×   | *   |
| D.02)               | 0  | 0   | 0   | 0   | 0   | *   | ×   | *   | *   | *   | ×   | *   |
| E.0                 | X  | *   | 0   | 0   | 0   | 0   | *   | ×   | X   | *   | ×   | ×   |
| F.0                 | X  | *   | 0   | 0   | 0   | 0   | ×   | ×   | ×   | *   | ×   | *   |
| G.0                 | ×  | *   | ×   | X   | ×   | ×   | 0   | ×   | ×   | *   | ×   | ×.  |
| H.x                 | ×  | ×   | ×   | *   | ×   | *   | *   | 0   | ×   | **  | ×   | *   |
| J.x                 | ×  | *   | ×   | ×   | ×   | *   | ×   | ×   | 0   | *   | ×   | ж   |
| K.2                 | ×  | *   | ×   | ×   | ×   | *   | ×   | *   | ×   | 0   | ×   | *   |
| L.0                 | ×  | ×   | X   | ×   | ×   | ×   | ×   | ×   | ×   | *   | 0   |     |
| M.0                 | ×  | *   | ×   | ×   | ×   | ×   | ×   | *   | ×   | *   | *   | 0   |

2) It is recommended that any FMS connected to an Rel D.0 or higher MP90 system should have at least revision C.0 or higher